

TIMM17B

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:18kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TIMM17B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

JM3; TIM17B; DXS9822;

Catalog #:A8477

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10245

Isotype:IgG

Swiss Prot:O60830

Purity:Affinity purification

For research use only.

Background:

This gene encodes a multipass transmembrane protein that forms an integral component of the mitochondrial translocase TIM23 complex. This complex facilitates the transport of mitochondrial proteins from the cytosol across the mitochondrial inner membrane and into the mitochondrion. There is a pseudogene for this gene on chromosome 12. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).