

NEK6

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:20 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NEK6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SID6-1512;

Catalog #: A8481

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10783

Isotype: IgG

Swiss Prot: Q9HC98

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a kinase required for progression through the metaphase portion of mitosis. Inhibition of the encoded protein can lead to apoptosis. This protein also can enhance tumorigenesis by suppressing tumor cell senescence. Several transcript variants encoding a few different isoforms have been found for this gene.

To place an order, please [Click HERE](#).