

SYNPO

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:99kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human SYNPO

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A8484

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11346

Isotype:IgG

Swiss Prot:Q8N3V7

Purity:Affinity purification

For research use only.

Background:

Synaptopodin is an actin-associated protein that may play a role in actin-based cell shape and motility. The name synaptopodin derives from the protein's associations with postsynaptic densities and dendritic spines and with renal podocytes (Mundel et al., 1997)

To place an order, please [Click HERE](#).