

LPIN1

Reactivity:Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:99kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LPIN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAP1;

Catalog #:A8486

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23175

Isotype:IgG

Swiss Prot:Q14693

Purity:Affinity purification

For research use only.

Background:

This gene encodes a magnesium-ion-dependent phosphatidic acid phosphohydrolase enzyme that catalyzes the penultimate step in triglyceride synthesis including the dephosphorylation of phosphatidic acid to yield diacylglycerol. Expression of this gene is required for adipocyte differentiation and it also functions as a nuclear transcriptional coactivator with some peroxisome proliferator-activated receptors to modulate expression of other genes involved in lipid metabolism. Mutations in this gene are associated with metabolic syndrome, type 2 diabetes, and autosomal recessive acute recurrent myoglobinuria (ARARM). This gene is also a candidate for several human lipodystrophy syndromes. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional splice variants have been described but their full-length structures have not been determined.

To place an order, please [Click HERE](#).