

TIMM10B

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 12kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TIMM10B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FXC1; Tim9b; TIM10B;

Catalog #: A8488

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26515

Isotype: IgG

Swiss Prot: Q9Y5J6

Purity: Affinity purification

For research use only.

Background:

FXC1, or TIMM10B, belongs to a family of evolutionarily conserved proteins that are organized in heterooligomeric complexes in the mitochondrial intermembrane space. These proteins mediate the import and insertion of hydrophobic membrane proteins into the mitochondrial inner membrane.

To place an order, please [Click HERE](#).