

TCEB3B

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 84kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TCEB3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ELOA2; HsT832; TCEB3L;

Catalog #: A8492

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51224

Isotype: IgG

Swiss Prot: Q8IYF1

Purity: Affinity purification

For research use only.

Background:

This gene encodes the transcriptionally active subunit of the SIII (or elongin) transcription elongation factor complex, which also includes two regulatory subunits, elongins B and C. This complex acts to increase the rate of RNA chain elongation by RNA polymerase II by suppressing transient pausing of the polymerase at many sites along the DNA template. Whereas a related protein with similar function, elongin A, is ubiquitously expressed, the encoded protein is specifically expressed in the testis, suggesting it may have a role in spermatogenesis.

To place an order, please [Click HERE](#).