

EMCN

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EMCN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EMCN2; MUC14;

Catalog #: A8494

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51705

Isotype: IgG

Swiss Prot: Q9ULC0

Purity: Affinity purification

For research use only.

Background:

EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al., 2001

To place an order, please [Click HERE](#).