

FAM20A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FAM20A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AI1G; AIGFS; FP2747;

Catalog #: A8496

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54757

Isotype: IgG

Swiss Prot: Q96MK3

Purity: Affinity purification

For research use only.

Background:

This locus encodes a protein that is likely secreted and may function in hematopoiesis. A mutation at this locus has been associated with amelogenesis imperfecta and gingival hyperplasia syndrome. Alternatively spliced transcript variants have been identified.

To place an order, please [Click HERE](#).