

SMEK1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 95kDa

Observed MW: Refer to figures

Immunogen:

Recombinant Protein of human SMEK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

smk1; FLFL1; PP4R3; SMEK1; smk-1; PP4R3A; MSTP033; KIAA2010;

To place an order, please [Click HERE](#).

Catalog #: A8500

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55671

Isotype: IgG

Swiss Prot: Q6IN85

Purity: Affinity purification

For research use only.