

RARS2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 66kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RARS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PCH6; ArgRS; RARSL; DALRD2; PRO1992; dJ382I10.6;

Catalog #: A8503

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57038

Isotype: IgG

Swiss Prot: Q5T160

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an arginyl-tRNA synthetase that is found in the mitochondrial matrix. Defects in this gene are a cause of pontocerebellar hypoplasia type 6 (PCH6).

To place an order, please [Click HERE](#).