

PPP4R4

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:99kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PPP4R4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PP4R4; CFAP14; KIAA1622;

Catalog #:A8505

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57718

Isotype:IgG

Swiss Prot:Q6NUP7

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a HEAT-like repeat-containing protein. The HEAT repeat is a tandemly repeated, 37-47 amino acid long module occurring in a number of cytoplasmic proteins. Arrays of HEAT repeats form a rod-like helical structure and appear to function as protein-protein interaction surfaces. The repeat-containing region of this protein has some similarity to the constant regulatory domain of the protein phosphatase 2A PR65/A subunit. The function of this particular gene product has not been determined. Alternative splicing has been observed for this gene and two transcript variants encoding distinct isoforms have been identified.

To place an order, please [Click HERE](#).