

TMEM43

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TMEM43

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LUMA; ARVC5; ARVD5; EDM7;

Catalog #: A8509

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79188

Isotype: IgG

Swiss Prot: Q9BTV4

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the TMEM43 family. Defects in this gene are the cause of familial arrhythmogenic right ventricular dysplasia type 5 (ARVD5), also known as arrhythmogenic right ventricular cardiomyopathy type 5 (ARVC5). Arrhythmogenic right ventricular dysplasia is an inherited disorder, often involving both ventricles, and is characterized by ventricular tachycardia, heart failure, sudden cardiac death, and fibrofatty replacement of cardiomyocytes. This gene contains a response element for PPAR gamma (an adipogenic transcription factor), which may explain the fibrofatty replacement of the myocardium, a characteristic pathological finding in ARVC.

To place an order, please [Click HERE](#).