

TEFM

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:42kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TEFM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C17orf42;

To place an order, please [Click HERE](#).

Catalog #:A8511

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79736

Isotype:IgG

Swiss Prot:Q96QE5

Purity:Affinity purification

For research use only.