

NUDT18

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 36kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NUDT18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MTH3;

Catalog #: A8514

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79873

Isotype: IgG

Swiss Prot: Q6ZVK8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the Nudix hydrolase family. Nudix hydrolases eliminate potentially toxic nucleotide metabolites from the cell and regulate the concentrations and availability of many different nucleotide substrates, cofactors, and signaling molecules. This protein contains a Nudix hydrolase domain and hydrolyzes oxidized forms of guanosine and deoxyguanosine diphosphates.

To place an order, please [Click HERE](#).