

CNTNAP3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 141kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CNTNAP3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CASPR3; CNTNAP3A;

Catalog #: A8515

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79937

Isotype: IgG

Swiss Prot: E9PY62

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the NCP family of cell-recognition molecules. This family represents a distinct subgroup of the neurexins. NCP proteins mediate neuron-glial interactions in vertebrates and glial-glial contact in invertebrates. The protein encoded by this gene may play a role in cell recognition within the nervous system. Alternatively spliced transcript variants encoding different isoforms have been described but their biological nature has not been determined.

To place an order, please [Click HERE](#).