

NAA25

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 112kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NAA25

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NAP1; MDM20; C12orf30;

Catalog #: A8516

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80018

Isotype: IgG

Swiss Prot: Q14CX7

Purity: Affinity purification

For research use only.

Background:

This gene encodes the auxiliary subunit of the heteromeric N-terminal acetyltransferase B complex. This complex acetylates methionine residues that are followed by acidic or asparagine residues.

To place an order, please [Click HERE](#).