

RNF34

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 42kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RNF34

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RFI; RIF; RIFF; hRFI; CARP1; CARP-1;

Catalog #: A8517

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80196

Isotype: IgG

Swiss Prot: Q969K3

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a RINF finger, a motif known to be involved in protein-protein and protein-DNA interactions. This protein interacts with DNAJA3/hTid-1, which is a DnaJ protein reported to function as a modulator of apoptosis. Overexpression of this gene in HeLa cells was shown to confer the resistance to TNF-alpha induced apoptosis, suggesting an anti-apoptotic function of this protein. This protein can be cleaved by caspase-3 during the induction of apoptosis. This protein also targets p53 and phospho-p53 for degradation. Alternatively splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).