

ANKRD44

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 108kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ANKRD44

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PP6-ARS-B;

To place an order, please [Click HERE](#).

Catalog #: A8521

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 91526

Isotype: IgG

Swiss Prot: Q8N8A2

Purity: Affinity purification

For research use only.