

GPRC6A

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 105kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GPRC6A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GPCR; bA86F4.3;

Catalog #: A8525

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 222545

Isotype: IgG

Swiss Prot: Q5T6X5

Purity: Affinity purification

For research use only.

Background:

Members of family C of the G protein-coupled receptor (GPCR) superfamily, such as GPRC6A, are characterized by an evolutionarily conserved amino acid-sensing motif linked to an intramembranous 7-transmembrane loop region. Several members of GPCR family C, including GPRC6A, also have a long N-terminal domain (summary by Pi et al., 2005)

To place an order, please [Click HERE](#).