

ULK1

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:1000 IF 1:10 - 1:100

Calculated MW: 113kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ULK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KIAA0722;

To place an order, please [Click HERE](#).

Catalog #: A8529

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8408

Isotype: IgG

Swiss Prot: O75385

Purity: Affinity purification

For research use only.