

LAMP2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:20 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human LAMP2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8531

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 3920

Isotype: IgG

Swiss Prot: P13473

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins.

To place an order, please [Click HERE](#).