

AQP8

Reactivity: Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human AQP8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8539

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 343

Isotype: IgG

Swiss Prot: O94778

Purity: Affinity purification

For research use only.

Background:

Aquaporin 8 (AQP8) is a water channel protein. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). Aquaporin 8 mRNA is found in pancreas and colon but not other tissues. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).