

GLP1R

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GLP1R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8547

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2740

Isotype: IgG

Swiss Prot: P43220

Purity: Affinity purification

For research use only.

Background:

GLP1R (glucagon-like peptide 1 receptor) is a protein-coding gene. Diseases associated with GLP1R include fasting hypoglycemia, and type 2 diabetes mellitus. GO annotations related to this gene include transmembrane signaling receptor activity and glucagon receptor activity. An important paralog of this gene is GHRHR.

To place an order, please [Click HERE](#).