

MYO5B

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 214kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MYO5B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8558

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4645

Isotype: IgG

Swiss Prot: Q9ULV0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene, together with other proteins, may be involved in plasma membrane recycling. Mutations in this gene are associated with microvillous inclusion disease.

To place an order, please [Click HERE](#).