

BCAN

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 99kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human BCAN

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9368

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 63827

Isotype: IgG

Swiss Prot: Q96GW7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]

To place an order, please [Click HERE](#).