

NUP107

Reactivity: Human Mouse Rat

Tested applications: WB IF IP

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 106kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human NUP107

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9373

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57122

Isotype: IgG

Swiss Prot: P57740

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the nucleoporin family. The protein is localized to the nuclear rim and is an essential component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).