

PYGM

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 97kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PYGM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9392

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5837

Isotype: IgG

Swiss Prot: P11217

Purity: Affinity purification

For research use only.

Background:

This gene encodes a muscle enzyme involved in glycogenolysis. Highly similar enzymes encoded by different genes are found in liver and brain. Mutations in this gene are associated with McArdle disease (myophosphorylase deficiency), a glycogen storage disease of muscle. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).