

THRAP3

Reactivity: Human

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 109kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human THRAP3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A9396

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9967

Isotype: IgG

Swiss Prot: Q9Y2W1

Purity: Affinity purification

For research use only.