

NUP35

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human NUP35

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9398

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 129401

Isotype: IgG

Swiss Prot: Q8NFH5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the nucleoporin family. The encoded protein contains two membrane binding regions, is localized to the nuclear rim, and is part of the nuclear pore complex. All molecules entering or leaving the nucleus either diffuse through or are actively transported by the nuclear pore complex. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene have been defined on chromosomes 7 and 10. [provided by RefSeq, Dec 2013]

To place an order, please [Click HERE](#).