

TMEM70

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 29kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TMEM70

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9423

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54968

Isotype: IgG

Swiss Prot: Q9BUB7

Purity: Affinity purification

For research use only.

Background:

This gene likely encodes a mitochondrial membrane protein. The encoded protein may play a role in biogenesis of mitochondrial ATP synthase. Mutations in this gene have been associated with neonatal mitochondrial encephalomyopathy due to ATP synthase deficiency. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2010]

To place an order, please [Click HERE](#).