

EFS

Reactivity: Human Mouse Rat

Tested applications: WB IF IP

Recommended Dilution: WB 1:500 - 1:2000 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human EFS

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9429

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10278

Isotype: IgG

Swiss Prot: O43281

Purity: Affinity purification

For research use only.

Background:

The longest protein isoform encoded by this gene contains an SH3 domain, which is known to be important in intracellular signal transduction. The protein encoded by a similar gene in mice was shown to bind to SH3 domains of protein-tyrosine kinases. The function of this gene is unknown. Three alternatively spliced variants encoding different isoforms have been described for this gene. [provided by RefSeq, Mar 2013]

To place an order, please [Click HERE](#).