

TNPO1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 102kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TNPO1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9435

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3842

Isotype: IgG

Swiss Prot: Q92973

Purity: Affinity purification

For research use only.

Background:

This gene encodes the beta subunit of the karyopherin receptor complex which interacts with nuclear localization signals to target nuclear proteins to the nucleus. The karyopherin receptor complex is a heterodimer of an alpha subunit which recognizes the nuclear localization signal and a beta subunit which docks the complex at nucleoporins. Alternate splicing of this gene results in two transcript variants encoding different proteins. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).