

T

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human T

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TFT; SAVA

Catalog #: A9439

Antibody Type:

Monoclonal Antibody

Species: 23

Gene ID: 6862

Isotype: IgG

Swiss Prot: O15178

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an embryonic nuclear transcription factor that binds to a specific DNA element, the palindromic T-site. It binds through a region in its N-terminus, called the T-box, and effects transcription of genes required for mesoderm formation and differentiation. The protein is localized to notochord-derived cells. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).