

MAZ

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:200 - 1:2000 IF 1:50 - 1:100

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human MAZ

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A9451

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4150

Isotype: IgG

Swiss Prot: P56270

Purity: Affinity purification

For research use only.