

CHCHD4

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100 IP 1:20 - 1:50

Calculated MW: 16kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human CHCHD4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9453

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 131474

Isotype: IgG

Swiss Prot: Q8N4Q1

Purity: Affinity purification

For research use only.

Background:

CHCHD4, a component of human mitochondria, belongs to a protein family whose members share 6 highly conserved cysteine residues constituting a -CXC-CX(9)C-CX(9)C- motif in the C terminus (Hofmann et al., 2005 [PubMed 16185709]). [supplied by OMIM, Mar 2008]

To place an order, please [Click HERE](#).