

ACE2

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 92kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human ACE2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9454

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 59272

Isotype: IgG

Swiss Prot: Q9BYF1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the angiotensin-converting enzyme family of dipeptidyl carboxydipeptidases and has considerable homology to human angiotensin 1 converting enzyme. This secreted protein catalyzes the cleavage of angiotensin I into angiotensin 1-9, and angiotensin II into the vasodilator angiotensin 1-7. The organ- and cell-specific expression of this gene suggests that it may play a role in the regulation of cardiovascular and renal function, as well as fertility. In addition, the encoded protein is a functional receptor for the spike glycoprotein of the human coronaviruses SARS and HCoV-NL63. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).