

SLC34A2

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC34A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9460

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10568

Isotype: IgG

Swiss Prot: O95436

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a pH-sensitive sodium-dependent phosphate transporter. Phosphate uptake is increased at lower pH. Defects in this gene are a cause of pulmonary alveolar microlithiasis. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, May 2010]

To place an order, please [Click HERE](#).