

FZD5

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 65kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human FZD5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9469

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7855

Isotype: IgG

Swiss Prot: Q13467

Purity: Affinity purification

For research use only.

Background:

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD5 protein is believed to be the receptor for the Wnt5A ligand. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).