

DLL4

Reactivity: Human Mouse

Tested applications: WB IP

Recommended Dilution: WB 1:200 - 1:1000 IP 1:20 - 1:50

Calculated MW: 75kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human DLL4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9471

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54567

Isotype: IgG

Swiss Prot: Q9NR61

Purity: Affinity purification

For research use only.

Background:

This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.
[provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).