

TCERG1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 124kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TCERG1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9479

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10915

Isotype: IgG

Swiss Prot: O14776

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear protein that regulates transcriptional elongation and pre-mRNA splicing. The encoded protein interacts with the hyperphosphorylated C-terminal domain of RNA polymerase II via multiple FF domains, and with the pre-mRNA splicing factor SF1 via a WW domain. Alternative splicing results in multiple transcripts variants encoding different isoforms. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).