

AHI1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 137kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human AHI1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9493

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54806

Isotype: IgG

Swiss Prot: Q8N157

Purity: Affinity purification

For research use only.

Background:

This gene is apparently required for both cerebellar and cortical development in humans. This gene mutations cause specific forms of Joubert syndrome-related disorders. Joubert syndrome (JS) is a recessively inherited developmental brain disorder with several identified causative chromosomal loci. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Oct 2008]

To place an order, please [Click HERE](#).