

ISG20

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human ISG20

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD25; HEM45

Catalog #: A9498

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3669

Isotype: IgG

Swiss Prot: Q96AZ6

Purity: Affinity purification

For research use only.

Background:

ISG20 (interferon stimulated exonuclease gene 20kDa) is a protein-coding gene. Diseases associated with ISG20 include malignant glioma, and leukemia. GO annotations related to this gene include single-stranded DNA 3'-5' exodeoxyribonuclease activity and exonuclease activity. An important paralog of this gene is REXO4.

To place an order, please [Click HERE](#).