

HES5

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:20 - 1:200

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human HES5

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9535

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 388585

Isotype: IgG

Swiss Prot: Q5TA89

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of basic helix-loop-helix transcriptional repressors. The protein product of this gene, which is activated downstream of the Notch pathway, regulates cell differentiation in multiple tissues. Disruptions in the normal expression of this gene have been associated with developmental diseases and cancer. [provided by RefSeq, Dec 2008]

To place an order, please [Click HERE](#).