

PHLPP1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 185kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PHLPP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCOP; PHLPP; PLEKHE1

Catalog #: A9542

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23239

Isotype: IgG

Swiss Prot: O60346

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the serine/threonine phosphatase family. The encoded protein promotes apoptosis by dephosphorylating and inactivating the serine/threonine kinase Akt, and functions as a tumor suppressor in multiple types of cancer. Increased expression of this gene may also play a role in obesity and type 2 diabetes by interfering with Akt-mediated insulin signaling.

To place an order, please [Click HERE](#).