

OXTR

Reactivity: Human

Tested applications: WB IHC IP FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human OXTR

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9550

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5021

Isotype: IgG

Swiss Prot: P30559

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the G-protein coupled receptor family and acts as a receptor for oxytocin. Its activity is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system. The oxytocin-oxytocin receptor system plays an important role in the uterus during parturition. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).