

EIF3K

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EIF3K

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9969

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27335

Isotype: IgG

Swiss Prot: Q9UBQ5

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).