

TRIM39

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TRIM39

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9977

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56658

Isotype: IgG

Swiss Prot: Q9HCM9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The function of this protein has not been identified. This gene lies within the major histocompatibility complex class I region on chromosome 6. Alternate splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).