

SAV1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SAV1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SAV; WW45; WWP4

Catalog #: A9980

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 60485

Isotype: IgG

Swiss Prot: Q9H4B6

Purity: Affinity purification

For research use only.

Background:

WW domain-containing proteins are found in all eukaryotes and play an important role in the regulation of a wide variety of cellular functions such as protein degradation, transcription, and RNA splicing. This gene encodes a protein with two WW domains, a SARAH domain, and a coiled-coil region and is ubiquitously expressed in adult tissues. This protein binds to MST1 (mammalian sterile 20-like kinase 1) and promotes MST1-induced apoptosis. It has also been shown to bind to HAX1 (hematopoietic cell-specific protein 1 (HS1)-associated protein X-1) and to attenuate the anti-apoptotic effects of HAX1. Studies in human and mouse suggest this gene acts as a tumor suppressor.

To place an order, please [Click HERE](#).