

UQCC2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UQCC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9986

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84300

Isotype: IgG

Swiss Prot: Q9BRT2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nucleoid protein localized to the mitochondria inner membrane. The encoded protein affects regulation of insulin secretion, mitochondrial ATP production, and myogenesis through modulation of mitochondrial respiratory chain activity. [provided by RefSeq, Oct 2012]

To place an order, please [Click HERE](#).