

TIMM50

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TIMM50

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A9991

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92609

Isotype: IgG

Swiss Prot: Q3ZCQ8

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).